




**NEW TECHNOLOGIES:
BIOMATERIALS, BIOMECHANICS,
BIOPRINTING**

Lisa Ferrara, Ph.D.



OrthoKinetic Technologies, LLC
&
OrthoKinetic Testing Technologies, LLC



DISCLOSURES

- OKT – Strategic Planning & Regulatory Consulting
- OKT² – ISO 17025 A2LA Accredited Test Facility



OVERVIEW

- What's new in biomaterials
- Mechanotransduction
- Combining Biomaterials, Nanostructures & Surfaces
- New & Exciting Technology



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

BIOMATERIALS



- A material intended to interface with biological systems to augment, treat, or replace damaged tissue

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

What's new in Biomaterials



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

PEEK

- Conventional
- HA PEEK
- HA Coated
- Titanium Coated
 - plasma spray, molded, printed inserts
- Coated
 - Many different Configurations
 - Endplates or Total Coverage
 - Printed, machined endplates, inserts
- Porous
 - Grown
 - Printed
 - Formed



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

POROUS PEEK

- DUKE & Georgia Tech.
- PEEK scoria grown directly from the bulk PEEK on surface
- Surface treatment Vertera Spine's COHERE system.

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

HA PEEK – InVibio – 10 Weeks Post-OP

Talos®-C (HA) IBF
2 Level HA PEEK and 2 Level Standard PEEK

Courtesy of Meditech Advisors for the Talos-C

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

HA COATED PEEK (BIOMIMETIC)

- Short half life of HA, Debonding issues
- Increase surface roughness on PEEK, CFR-PEEK – biomimetic process
- Simulated Body Fluids / CO2 supply vs. Plasma Sprayed HA
- Increased Alk Phos. Activity
- Increased osteogenic differentiation

Figure 2. ALP expression normalized by the amount of DNA in each sample at day 7, 14 and 21.


Material	Day 7	Day 14	Day 21
PEEK	~0.010	~0.015	~0.020
CFR-PEEK	~0.015	~0.020	~0.025
PLGA	~0.010	~0.015	~0.020
PLGA-coated HA	~0.020	~0.030	~0.040

A biomimetic hydroxyapatite coating on polyetheretherketone
Sara Amini, Gordon Blum, and Melanie Coakley
*University College London, Institute of Orthopaedics and Musculo-Skeletal Science, United Kingdom

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

WHAT IS ADDITIVE MANUFACTURING?

- Process of joining materials to make objects from 3D models
- Not subtractive – AM adds material in layers
- Builds layer upon layer onto substrates
- Value added
 - Can build complex structures that are difficult to machine
 - Characterized by rough surface quality – improved cellular response
 - Designed, controlled, interconnected porosity
 - High purity – restricted presence of O₂
 - High power processing velocity – rapid serial production of implants
 - Customize fit to patient anatomy



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

3D PRINTING - ADDITIVE MANUFACTURING

- Vehicle to fabricate open, complex structures for improved interface responses
- Microstructures, working on perfecting nanosizes
- Design custom porosities - graded porosities & structures
- Build Nanostructures & Nanoscaffolds
- Different material = different mechanical stimuli & response
- Programmed responses (mechanotransduction at cell level)



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

Additive Manufactured Implants


- Biomaterial + Manufactured Process
 - Titanium Structures
 - CoCr
 - PEEK
 - Ceramics
- Provides a VEHICLE for fabricating complex devices, osseointegrative surfaces, graded materials w/o bonding processes



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

PEEK & CAD CAM

- Additive Manufacturing
- Early Stage
- PEEK partial denture frameworks
- Strong & lightweight improved patient comfort;
- Digital design matches the patient's individual anatomy
- Taste-neutral
- Fabricate restorations



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**


Titanium: Printed Coated or AM

- AM or 3D Printing is not the holy grail
- Mechanical behavior that is the **key**
 - Stress transfer
 - Strain behavior
 - LOCAL MODULUS vs. BULK MODULUS
- Failure Mode - IMPORTANT




ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

Cspine - Sagittal View




- No evidence of graft material at the levels with PEEK cage (yellow arrows)
- Partial consolidation and graft maturation at 4-Web level (blue arrows)
- 1 year out
- Courtesy of Dr. Camissa & 4Web



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

Multi-Level – Cspine / LSpine

3-Level



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

COMBINING BIOMATERIALS WITH INNOVATIVE DESIGN & SURFACE TREATMENTS



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

IMPROVING INTERFACE MECHANICS

- Why?
- Balance between Time & Stabilization -Earlier Stabilization = Earlier Fusion
- Increase time to heal = increase risk of failure
- Earlier stabilization = Decreased time to heal & Failure risk
- Improvements in materials & surfaces
- Improved Interface mechanics



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

INTERFACE MECHANICS

- Improve designs - Multi-directional bone exchange = better biomechanical stability
 - Full open, multiplanar
 - Large windows, multiple directions
 - Surface textures – endplates or entire device
 - Materials that promote osseointegration, osteoinduction, osteoconduction
- Improved stability allows for further ingrowth towards center- throughout implant
- Greater area of bone contact, greater fusion footprint – better stress distribution

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

Programmable Structures: Scaffolds, Trusses & Surface Modulation

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**


IMPORTANT VARIABLES


- Delivery
 - Cell Suspensions
 - Tissue-like constructs (scaffolds)
- Chemical properties
 - Growth factors
 - Degradation particles
 - ECM surface
- Physical properties
 - Structure
 - Topography
 - Rigidity
 - Mechanical Loading

Modify Cell Behavior

Survival
Organization
Migration
Proliferation
Differentiation

Optimize Cellular Response





Can Tissue Regeneration be 'programmed'?


GEOMETRICAL FACTORS-
OPTIMAL POROSITIES FOR BONE GROWTH


- Bone tissue penetration – pore sizes 200-900µm (ideal 300-500µm)
- Enhanced cell penetration through pores
- Enhanced extracellular matrix production


- Vertical pores 300µm
- Horizontal pores 200-400µm
- Surface Posts - 10µm x 10µm / 10µm spacing
 - Kim et al. - Enhanced growth for CTPs

}

3D Scaffold + Surface
Microtexture











Surface Textures - Mechanotransduction

- Induce mechanotransduction
- Guide cellular growth
- Increased cell proliferation – or -
- Increased cell efficiency
- On Metals: use of surface textures
 - Increase fatigue performance of implant
 - Careful not to compromise strength
 - Titanium = improve osseointerdigitation











MECHANOTRANSDUCTION PROCESS


- Process where cells sense & respond to mechanical stimuli
- Cells convert energy to electrical & biochemical signals that elicit response.
- Stem cells sense mechanical loading thru membrane-anchored mechanosensors
- Converted to biochemical signals
- Trigger multi-step activation in array of signaling cascades in cytoplasm
- Activation of transcription factor
- Translocate to nucleus and modulate expression of mechanosensitive genes










OVERVIEW






- Cells sense mechanical forces -convert into biological responses by triggering restructure of cytoskeleton
- FEEDBACK SYSTEM
- Cells respond







EXAMPLE OF MECHANOTRANSDUCTION

- Ti Coated
- Centinel Spine - ACTIVE™
- Cell migration of cultured cells studied by scanning electron microscopy (SEM).
- PEEK- absence of any cell attachment to surface.
- Ti-ACTIVE™ - coated device shows cell attachment
- Ti-ACTIVE™ - cell migration seen by cytoplasmic expansions of cultured cells (purple)




<http://www.centinelspine.com>





DIFFERENTIATION MSCs INTO SPECIFIC CELL PHENOTYPES

- Mechanical Forces
 - Varied/controlled stiffness of substrate
 - Differentiate with different stiffnesses
- MSCs respond to hydrogenated amorphous carbon nanotopographics with groove or grid surfaces
 - Groove structures demonstrated dynamic effect assoc. w/ stem cell alignment + elongation.
- Surface topography – Nanogrooves
 - Width/spacing – 40/30µm – optimal
 - Induced hMSCs -acquire neuronal characteristics in absence of differentiation agents.
- Mechanics = differentiation




Mechanotransduction: use the force(s)
 Paluch et al. BMC Biology (2015) 13:47
 DOI 10.1186/s12915-015-0150-4

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

SURFACE STRAIN


- Surface roughness had micro & nano structures that form peaks / valleys - strain responses during loading – micro and nanostrains measured
- Study showed that optimal biological responses to specific surface geometries Lending to different strain patterns



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**


SCAFFOLD STRUCTURES - VARIABLES AFFECT CELL RESPONSES

- Surface chemistry
- Matrix topography
 - Cell organization, alignment
 - Fiber alignment -> tissue development
- Rigidity
- Porosity
 - Large interconnected
 - small disconnected



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

MICRO-NANO SCAFFOLDS CELL DIFFERENTIATION



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**


MICROTEXTURED POLYMERS

- MEMS structures – soft lithographic process
- Assessed osteogenesis- proliferation, differentiation
- Connective tissue progenitor cells [CTPs] – Human Bone Marrow
- Cultured 60 days
 - Smooth PDMS
 - Microtextured w/Posts PDMS - 10µm dia/ 6µm height / 10µm spacing
- DNA quantification
- Western Blot Analysis – for Integrin α5 expression
- Real time RTP chain rxn – alk phosphatase gene expression

Post microtextures accelerate cell proliferation and osteogenesis

Eun Jung Kim^{1,2}, Cynthia A. Boehm^{1,2}, Alvaro Mata¹, Aaron J. Fleischman¹, George F. Muschler^{1,2}, Shuvo Roy^{1,2*}

Acta Biomaterialia 6 (2010) 160-169




ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

TEXTURED SURFACES -TISSUE SCAFFOLDS

- Disc & ligament regeneration
- Organized / programmed growth
- Respond to mechanical, chemical, electrical stimuli


Directed growth of a fibroblast that responds to mechanical stimuli



Astrocyte cells grown on
 (a) non-woven PET fiber matrix
 (b) microfabricated PLGA scaffold


ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

RESULTS





	Control (Smooth PDMS)	Post Microtextures
Day 9	0.298	1.100
Day 30	0.347	1.094
Day 60	1.021	1.846

- Gene Expression for Collagen Type I




- Cell Proliferation & Growth Orientation





BIOACTIVE SURFACE TOPOGRAPHIES

- 3D Microstructure – defined microarchitecture + microtextures
- Provide specific physical cues to cells / tissue
- Precisely engineered architecture
- Tailored surface topography
 - Optimized for cell type
 - Microfabrication + Soft Lithography to build scaffolds
- Connective tissue progenitor cells [CTPs] – 3 Patients
- Cultured **9** days under osteoblast phenotype condition
 - Smooth PDMS 3D structure
 - Microtextured w/Posts PDMS - 10µmX10µm
 - 66% Porosity
 - 5 layered scaffolds – 3D





200 µm
Textured




200 µm
Smooth

Biomaterials 3D (2009) 4610-4617
A three-dimensional scaffold with precise micro-architecture and surface micro-textures
Shawn Miao^{1,2,3,4}, Sun Jung Kim^{1,2}, Cynthia A. Barber^{1,2,3}, Aaron J. Heitschman^{1,2}, George R. Maniatty^{1,2,3,4}, Glenn Yang^{1,2,3,4}





3D SCAFFOLD + MICROPOSTS

- 5 layer PDMS 3D porous structure





800 µm




100 µm
200 µm

Microtexture →
Microposts →





RESULTS

Cell Number/Colony on 3D Scaffolds



Experiment	Texture	Smooth
experiment #1	~450	~180
experiment #2	~650	~200
experiment #3	~480	~200



FINDINGS - 3D TEXTURE & 3D SMOOTH

- Cells attached, proliferated, migrated vertically betw layers, differentiated
- All 3 dimensions on different features of scaffold
- Cells grew top down to 4 layers deep
- Cells grew top/ bottom of each layer, walls of pores, columns
- Vertical direction, cells proliferated and migrated up to 1.2mm
- Colonies reached 3mm dia. horizontally
- Cytoskeleton grew within & around posts

THE FUTURE

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**



SMART SCAFFOLDS

- Can create repeatable & predictable MICRO and NANO architectures for optimized cellular responses
- Can modulate, differentiate cells – surface textures and structures
- Can layer for smart scaffold – multilayer organized growth
- Mechanical behavior of scaffold / substrate material can dictate cellular response and organization (stiffness)
- Multiplanar stress transfer – optimized interface mechanics
- Create 'passively smart' structures for cell differentiation and tissue regeneration

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

WHAT IS ADDITIVE MANUFACTURING?


- Process of joining materials to make objects from 3D models
- Not subtractive – AM adds material in layers
- Builds layer upon layer onto substrates
- Value added
 - Can build complex structures that are difficult to machine
 - Characterized by rough surface quality – improved cellular response
 - Designed, controlled, interconnected porosity
 - High purity – restricted presence of O₂
 - High power processing velocity – rapid serial production of implants
 - Customize fit to patient anatomy



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

AM PROCESSES & BIOMIMICRY

- Can build structures & functions that mimic biology & complex shapes of tissues
 - Bone → Truss structures



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

BIOPRINTING




ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

BIOPRINTING

Morphing new technologies to obtain impactful innovation

Regenerative Medicine – Next Wave



3D Bioprinting + Nanotechnology + Stem Cells
>>> INNOVATION



ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

BIOPRINTING USING AM

- ORGANS / FACIAL PRINTING / TISSUE SCAFFOLDS
 - Artificially constructs living tissue / organs
 - Spheroids contain aggregate of tens of thousands of cells
 - Builds layer-upon-layer of living cells from the patient's own body
 - Reduces rejection to almost zero.
 - [Organovo](#) + [Invetech](#) created commercial bioprinter - [NovoGen MMX](#). Loaded with bioink spheroids w/cells, substrated collagen / protein matrix etc....
 - December 2010, [Organovo](#) created [first blood vessels bioprinted using cells cultured from a single human](#).




(A) Bioink spheroids printed into layer of bioappet gel
(B) Additional layers printed to build object
(C) Stack spheroids face together and bioappet discovers
(D) Final 3D printing tissue

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

3D PRINTED TISSUE

- 3-D printed CaP bone substitute
- 3D of Bone



The collage includes a 3D printer on the left, a grid of colored dots in the center, a blue arrow pointing right, and a 2x2 grid of images on the right: a 3D printed bone, a 3D printed mesh, a 3D printed pelvis, and a close-up of a bone substitute in a surgical setting.

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

3D PRINTED SKIN & CARTILAGE

- [Print skin directly onto burn victims](#)
- **BIOPRINTING** - computer-controlled machine assemble biological matter using organic inks and super-tough thermoplastics.
- HUMAN CARTILAGE - flexible mats of electrospun synthetic polymer combined, layer-by-layer, with cartilage cells from rabbit ear deposited by ink jet printer.
- The constructs were square with a 10cm diagonal and a 0.4mm thickness.
- Synthetic materials ensure the strength of the construct and natural gel materials provide an environment that promotes cell growth.




The collage includes a 3D printer on the left, a 3D printed skin construct in the center, and a 2x2 grid of images on the right: a 3D printed skin, a 3D printed cartilage, a 3D printed hand, and a 3D printed arm.

ORTHO KINETIC TECHNOLOGIES, LLC **ORTHO KINETIC TESTING TECHNOLOGIES, LLC**

3D PRINTED ORGANS

- [Scientists Have 3D-Printed Mini Human Livers for the First Time Ever](#)
- New mini livers – live 40 days
- Great approximations of real deal - producing the proteins that carry hormones & drugs throughout body, cholesterol, and the major detoxification enzymes that let you consume reasonable amounts of alcohol without dying.
- **Organovo** plans to move on to the normal-sized organs



The collage includes a 3D printer on the left, a 3D printed organ in the center, and a 2x2 grid of images on the right: a 3D printed organ, a 3D printed organ, a 3D printed hand, and a 3D printed arm.

OrthoKinetic Technologies, LLC **OrthoKinetic Testing Technologies, LLC**

Evolution of Technology & Medicine

- Continue to innovate
- Improve outcomes, quality of life, longevity
- Continue to progress – good & bad



OrthoKinetic Technologies, LLC **OrthoKinetic Testing Technologies, LLC**

THANK YOU



CONTACT INFORMATION:
Lisa Ferrara, Ph.D.
lisa@orthokintech.com
Voice: 910.253.9883
www.orthokintech.com

OrthoKinetic Technologies, LLC **OrthoKinetic Testing Technologies, LLC**

