

Spinous process/intra-laminar devices?

Ira Fedder, MD
Keywest, 2018


Spinous Process




Spinous Process



Figure 7.6. Posterior apposition from an unstable disc resulting in "kissing" spinous processes and possible bursa development.






Device purpose ?

- Spacer – passive/active
- Motion preserving device
- Fusion device- stand alone
- Fusion augmentation-ALIF,LLIF,TLIF





Inter-laminar vs spinous process



Deyo 2013

- Medicare data - 75 yo
- 2006-2009
- 6000 spacers vs 17000 fusions
- Reduce LOS (1.4 vs 2.7 days)
- Reduced complications (1.2% vs 3.3%)
- 2 yr revision (17% vs 10)

Biomechanical impact– Intact? Decompressed? Normal Spine?

- Reduce motion in extension
- Reduce intradiscal pressure in extension
- Reduce facet loads
- Prevents(?) loss of foraminal area in extension
- Other effects may be dependent on the size and placement of the implant

In vitro studies-Aspen

- stand alone decreases ROM in Ext/Flexion
- with ALIF equal to bilateral pedicle screws
- with ALIF superior to anterior plating
- created flexion at the treated level



Aspen/TLIF -- Kaibara et al

- Stand alone reduced flex/ext similar to bilat pedicle screws
- Stand alone inferior in lateral bending and rotation
- Aspen/TLIF and unilateral screws similar to TLIF and bilateral screws

Spinal Simplicity
Innovative simple solutions

Minuteman G3R

- Interspinous fusion device
- Bone graft window

Standard Minuteman G3 Implant

Spinal Simplicity
Innovative simple solutions

Clinical Benefits

- Single Position Anterior/Posterior Fusion
 - MIS Lateral surgical approach, no need to reposition patient
- Fusion Device
 - Promotes fusion by utilizing the large bone graft window and HA coated implants
- Immediate, Rigid Stabilization
 - Rigidity comparable to interbody fusion with supplemental pedicle screw-rod fixation
 - Supraspinous ligament remains intact acting as tension band

Spinal Simplicity
Innovative simple solutions

Biomechanics




Figure 1. Total relative motion between L3 and L4. All data are represented as mean \pm standard deviation.

- After LLIF, Minuteman provided:
 - similar stability to bilat. pedicle screws and facet screws in F-E and axial rotation

Spinal Simplicity
innovative simple solutions

Biomechanics

- The Minuteman interspinous fixation/fusion enhances stability of the lumbar spine similar to other posterior fixation techniques when used in a DLIF/XLIF construct
- The Minuteman device may be an effective minimally-invasive alternative to traditional pedicle screw/rod constructs when used to support a DLIF/XLIF

